

Preeclampsia risk stratification early in pregnancy:  
Conversion of a promising metabolomics discovery into a LC-MS based clinical assay

**Liz Bond**<sup>1,6</sup>, Charline Lenaerts<sup>2</sup>, Christopher Benton<sup>3</sup>, Phil Baker<sup>4,6</sup>, Louise Kenny<sup>5,6</sup>, Robin Tuytten<sup>1,6</sup>




<sup>1</sup>Metabolomic Diagnostics, Little Island, Cork, Ireland; <sup>2</sup>Lab. of Pharmaceutical Analysis, Faculty of Medicine & Pharmacy, University of Mons, Belgium; <sup>3</sup>Agilent Technologies, Life Sciences & Chemical Analysis Group, Stockport, UK; <sup>4</sup>Gravida: National Centre for Growth and Development, The University of Auckland, New Zealand; <sup>5</sup>The Irish Centre for Fetal and Neonatal Translational Research, Cork University Maternity Hospital, Cork, Ireland; <sup>6</sup>on behalf of the IMproved Pregnancy Outcomes by Early Detection consortium (IMPROVED)



## Abstract

Basic metabolomics research has uncovered that combinations of blood borne metabolites can risk-stratify women early in pregnancy according to their risk of developing preeclampsia later in their pregnancy. Since then, a company has been established which is dedicated to translating this finding into a tool for health care providers and pregnant women. A targeted approach is being developed whereby ca. 40 metabolites are (semi-) quantified using liquid chromatography-tandem mass spectrometry. An update on the method development progress as well as an overview of the clinical studies lined-up to verify and validate the preeclampsia risk stratification test will be discussed.



**Fig.1.** Overview Metabolomic Diagnostics' Preeclampsia Risk Stratification product.

## Background

- Preeclampsia is mostly a syndrome of late pregnancy characterised by concomitant (new onset of) hypertension (high blood pressure) accompanied by new onset proteinuria (elevated protein in the urine) or signs of multi-systemic impairment. Currently there is no cure for preeclampsia other than delivery.
- Preeclampsia has been the most significant cause of maternal death over recent decades: between 70,000 and 80,000 women die every year from preeclampsia and in excess of half a million new born infants die annually as a direct result of the condition.
- Efficient prediction of preeclampsia is considered a crucial step stone to deliver more effective prenatal care, minimize preeclampsia related complications and per result reduce health care costs.

## The Challenge


- Thus far preeclampsia prediction is largely depending on clinical risk factors, but these are marginally useful in healthy first time pregnant women who account for >50% of preeclampsia cases.
- Prediction of preeclampsia in first time pregnant women requires a panel of biomarkers in order to encapsulate the complex pathogenesis of the syndrome.
- Research and (commercial) test development requires prospectively collected 1<sup>st</sup> pregnancy biobanks, whereby 1000's of women need to be longitudinally monitored and sampled throughout pregnancy.

## The Opportunity

- Kenny *et al*<sup>1</sup> from the University College Cork, Ireland found that accurate prediction of preeclampsia in 1<sup>st</sup> time pregnant women is possible using a panel of blood-borne metabolite present in plasma of ~15 weeks pregnant women.
- In recent years there is an increased demand for the rationalisation of health-care: identification of 1<sup>st</sup> time women at risk of preeclampsia will allow health-care workers to better administer the right prenatal care to the right women.
- Metabolite analysis using LC-MS is well established in pharma (bioanalysis) and neonatal screening (“heel prick”). LC-MS is therefore an obvious platform choice to port a potentially disruptive preeclampsia risk stratification test into the clinical laboratory.
- Together the above present an entrepreneurial opportunity to develop a disruptive, first-in-class diagnostic tool to stratify 1<sup>st</sup> time pregnant women early in pregnancy according to their preeclampsia risk.

## The Approach

- Collaboration between Company, technology inventors and other clinical experts.
- Creation of Public-Private partnerships to lever transnational funding to propel both academic research and SME-based product development.
- Seminal partnership with dedicated instrument manufacturer.
- Understanding of funding mechanisms - what funding to apply / seek at what stage of company / product development.




**Fig.2.** Overview of the interplay between funding bodies, academia and company

## Method Development


### Selection of Target Metabolites

- Primarily based on metabolites identified in Kenny *et al.* (2010)<sup>1</sup>; 41 Target Metabolites selected
- Selection of corresponding stable isotopically labelled standards chosen based on the metabolite classes; 16 ISTDs selected



**Fig.3. Composition of metabolite panel**

## Development of multiplex LC-MRM assay for target metabolites and selection of ISTDs




LC:  
10 minute gradient elution -Agilent PFP Pursuit col;  
separation of hydrophilic and hydrophobic compounds

MS:  
Unique Qual and Quant transitions for 41 Target  
Metabolites and 16 ISTDs: 112 MRM transitions / run  
(dMRM)

### Development of a simple, comprehensive plasma extraction procedure


- Metabolite extraction procedure from plasma developed and optimised for all metabolites (and ISTDs) and compatible with LC-MS analysis methodology.



**Fig.4.** Metabolite sample processing pipeline

### Proof of Principle: Determination of target metabolites in pregnancy blood

- 41 metabolites of interest are assay-able in pregnancy blood (15wks) with good precision
- **Semi quantitative data** generated using Relative ISTD quantitation
- Current throughput: 435 samples/week (QC samples excluded)
- **Ongoing**: stability testing batch processing and reference range determinations



**Fig.5.** Precision metrics 40 metabolite assays; \* 5 independently prepared samples - %CV for whole sample processing pipeline

## The Near Future

- Case Control study - SCOPE (n ~700): Metabolite candidate confirmation and prediction algorithm refinement.
- Final Standard Operation Procedures for processing pipeline and data analysis.
- Field testing of SOPs and verification of algorithm in full SCOPE Europe (n ~2500).
- Technical and clinical validation of PrepSia v1.0 in IMPROVED (n ~5000)<sup>2,3</sup>.

## Conclusions

- Within Metabolomic Diagnostics, a single step metabolite extraction and a targeted LC-QqQ-MS approach using stable isotope labelled metabolites for relative quantification has been successfully developed.
- All major components are in place to commence processing the clinical samples as available in SCOPE and IMPROVED biobanks.

## References

1. Kenny, L. C. *et al.* Robust early pregnancy prediction of later preeclampsia using metabolomic biomarkers. *Hypertension* **56**, 741-9 (2010).
2. North, R. A. *et al.* Clinical risk prediction for preeclampsia in nulliparous women: development of model in international prospective cohort. *BMJ* **342**, d1875 (2011).  
([www.scopestudy.net](http://www.scopestudy.net))
3. Navaratnam, K. *et al.* A multi-centre phase IIa clinical study of predictive testing for preeclampsia: improved pregnancy outcomes via early detection (IMPROVED). *BMC Pregnancy Childbirth* **13**, 226 (2013).  
<http://www.fp7-improved.eu/>  
[www.metabolomicdiagnostics.com](http://www.metabolomicdiagnostics.com)

## Acknowledgements

The authors gratefully acknowledge funding from the EU-HEALTH Project IMPROVED (305169) of the Seventh Framework Programme (FP7).

Ongoing and insightful support from Agilent Technologies, Life Sciences & Chemical Analysis Group

